


Configuration Sheet

Commissions no./PatientID: _____


Step 1: What is ordered?

- evomove orthokit** (incl. holder for evomove®)
- evomove solokit** (incl. pocket)
- Shank
- left leg right leg
- Test supply (4 weeks)



Step 2 - 1st Option: 3D-Scan


3D-Scan by using a scanner.



Important: Please note when scanning, that knees and ankles are visible. Preferred .stl-file.

Step 2 - 2nd Option: Measuring with photos

U1	Knee (circumference)	cm
U2	Venter (circumference)	cm
U3	Tigh (circumference)	cm
L1	Lower leg (length)	cm



Three photos of the position of the electrodes and their size.

Important: The knees and ankles must be visible when taking pictures.

Step 3: Details of the Cuff

- Cable outlet mediale side laterale side


L3 Length
Top edge to cable outlet cm

U4 Distance
Tibial edge to the cable outlet cm

or

- Pocket medial lateral*

*Laterally only possible without tibial electrodes.



L2 Desired length
of the cuff cm

Notes:

Please send the completed configuration sheet with all pictures or scan-data to bestellung@evomotion.de.