


# Configuration Sheet Shank

Commissions no./PatientID: \_\_\_\_\_


## Step 1: What is ordered?

- evomove orthokit** (incl. holder for evomove®)
- evomove solokit** (incl. pocket)
- left leg  right leg
- Test supply (4 weeks)



## Step 2 - 1st Option: 3D-Scan


3D-Scan by using a scanner.



**Important:** Please note when scanning, that knees and ankles are visible. Preferred .stl-file.

## Step 2 - 2nd Option: Measuring with photos

U1	Knee (circumference)	cm
U2	Venter (circumference)	cm
U3	Tigh (circumference)	cm
L1	Lower leg (length)	cm



**Three photos** of the position of the electrodes and their size.

**Important:** The knees and ankles must be visible when taking pictures.

## Step 3: Details of the Cuff

- Cable outlet  mediale side  laterale side


L3 Length  
Top edge to cable outlet  cm

U4 Distance  
Tibial edge to the cable outlet  cm

or

- Pocket  medial  lateral\*

\*Laterally only possible without tibial electrodes.



L2 Desired length  
of the cuff  cm

Notes: \_\_\_\_\_

Please send the completed configuration sheet with all pictures or scan-data to [bestellung@evomotion.de](mailto:bestellung@evomotion.de).

# Configuration Sheet – Thigh


Commission No. / Patient ID:

## Step 1: What is ordered?

- evomove orthokit** (incl. holder)
- evomove solokit** (incl. pocket)
- Shorts  Cuff
- left leg  right leg
- female  male
- Test supply (4 weeks)


## Step 2 - 1st option: 3D scan

Create a **3D scan** using a scanner.




**Important:** Please note when scanning, that both legs can be seen standing from knee to waist. Preferably as an **.stl file**.

## Step 2 - 2nd option: Measurement with photos



**Important:**  
The seat should be straight and without upholstery.



Photos with appropriate electrodes.

**Important:**  
When taking photos, both legs must be visible from waist to knee.

**Do not forget:**  
Photograph from both sides laterally.

We recommend marking positions **U1**, **U2** and **U4** on the body for measuring. Measurements are made from mark to mark.

## Circumferences

<b>U1</b>	Waist (centered between the lower ribcage and the iliac crest, height of the belly button)	cm	
<b>U2</b>	Buttocks (strongest point buttocks/hip)	cm	
<b>Important:</b> Please measure both legs.		left	right
<b>U3</b>	Thigh (1cm below crotch)	cm	cm
<b>U4</b>	Knee circumference (back of the knee)	cm	cm

## Lengths


<b>L1</b>	Knee length (waist to knee circumference U4)	cm
<b>L2</b>	Buttocks (waist - strongest part of the buttocks)	cm
<b>L3</b>	Seat length (waist to seat)	cm

### Step 3: Details of the shorts/cuff


T1 Cell phone pocket left  right  none

T2 evomove® pocket lateral  medial

Notes:



T1



T2

evomove® pocket (lateral)

Shorts shown only schematically.

Please send the completed configuration sheet with all images or scan data to [bestellung@evomotion.de](mailto:bestellung@evomotion.de).