

Configuration Sheet – Shorts

Commission no. / Patient ID:
(kein vollständige Patient:innen Name!)


Reference (Ordernumber, patient ID ...):

Step 1: What is ordered?

evomove orthokit (incl. holder) evomove solokit (incl. pocket)	right leg left leg	female male	Fabric color (only one color possible) grey black	Test supply (4 Weeks) Continuous supply after test supply Continuous supply
Pair of new shorts (excl. Control unit)	Control unit (Optional) Hand (Smartphone) Arm (Watch)			

Schritt 2: Measurements (required in any case!)

Circumferences




U1 Waist (midway between the lower costal arch and the iliac crest, level with the navel)	cm
U2 Buttocks (strongest part buttocks/hips)	cm
U3 Thigh (1cm below crotch)	cm
U4 Knee circumference (back of the knee)	cm


Important:
Please measure both legs.

left

right

Lengths






Body dimensions are measured, e.g. using a rubber band to mark the positions U1, U2 and U4.

L1 Knee length (Waist to knee circumference U4)	cm
L2 Buttocks (Waist – strongest point on the buttocks)	cm
L3 Seat length (Waist to seat) Important: The seat should be straight and without cushions.	cm
Body height	cm

Step 3: Detection of the electrode position (Option, please check)

1. Option: 3D-Scan



Create a 3D scan using a scanner.
We accept .obj and .stl files.

Important: When scanning, please make sure that the buttocks and both legs can be seen standing from knee to waist.

2. Option: Photos





Photo with corresponding electrodes.
Important: When taking the photograph, both legs must be visible from waist to knee while standing. Make sure that the camera is aligned at the same height as the electrodes.

Do not forget: Photograph laterally from both sides.

Step 4: Details of the shorts

1	No			
Cord in waistband	Yes			
2	none			
Smartphone pocket	2a	right		
	2b	left		
3	No			
Zipper (Desired position may be adjusted depending on the electrode position)	Yes, in the waistband			
	3a	right	short	long
	3b	center	short	long
	3c	left	short	long
	Yes, on the leg			
	3d	right	short	long
	3e	left	short	long



T1 Pocket or cable outlet	Standard (lateral)	
	Individual (please mark on photo or scan)	
Wearing height	Waist shorts	Hipster shorts

Notes

Please send the **completed** configuration sheet with all images or scan data and any other comments to bestellung@evomotion.de.